

Export Chromatogram Consensus

To export consensus sequence use the *Consensus* tab of the *Options Panel*:

The screenshot shows the 'Consensus' tab of the 'Options Panel'. It is divided into two main sections:

- Consensus mode:**
 - Consensus type:** A dropdown menu currently set to 'Strict'.
 - Threshold:** A slider control set to 100%.
- Export consensus:** A section with a light red background containing:
 - Export to file:** A text input field with the value 'Aligned reads_consensus_3.txt' and a file selection icon (three dots).
 - File format:** A dropdown menu currently set to 'Plain text'.
 - Keep gaps:** An unchecked checkbox.
 - Export:** A button at the bottom of the section.

The following parameters are available:

Export to file - here you need to select path for the output file.

File format - format for the output file.

Keep gaps - set this parameter if the result consensus must keep gaps.

When you click on the *Export* button the consensus sequence will be exported into selected output file.