


Export Chromatogram Consensus

To export consensus sequence use the *Consensus* tab of the *Options Panel*:



The image shows a software interface for exporting a consensus sequence. It is divided into two main sections: 'Consensus mode' and 'Export consensus'. The 'Consensus mode' section includes a dropdown for 'Consensus type' set to 'Strict' and a 'Threshold' slider set to 100%. The 'Export consensus' section, highlighted in a light red background, includes a text input for 'Export to file' with the value 'Aligned reads_consensus_3.txt', a 'File format' dropdown set to 'Plain text', an unchecked checkbox for 'Keep gaps', and an 'Export' button.

The following parameters are available:

Export to file - here you need to select path for the output file.

File format - format for the output file.

Keep gaps - set this parameter if the result consensus must keep gaps.

When you click on the *Export* button the consensus sequence will be exported into selected output file.