

# Selecting Restriction Enzymes

To select an enzyme check it in the list. Notice that the enzyme appears in the *Selected enzymes* area of the dialog.

You can also use the *Select All* button to select all the enzymes available, the *Select None* button to deselect all the enzymes.

To select all enzymes with recognition sequence length shorter than the specified value click the *Select by length* button and input the minimum length in the dialog appeared.

To invert selection click the *Invert selection* button.

As soon as enzymes are selected you can click the *OK* button to search for corresponding restriction sites in the sequence.