


# Creating New Project

A *project* store links to the data files, cross-file data associations, and visualization settings.

Below is a description of how to create a new project manually. Note that if you have no project created when opening a file with a sequence, an alignment or any other biological data, a new anonymous project is created automatically.

To create a new project select the *File* *New project* menu or click the *New project* button on the main toolbar.

The dialog will appear:



The image shows a dialog box titled "Create New Project". It has two input fields. The first is labeled "Project name" and contains the text "New Project". The second is labeled "Save project to file" and contains the text "/home/..|project.uprj". To the right of the second field is a small button with three dots "...". At the bottom of the dialog, there are three buttons: "Help", "Cancel", and "Create".

Here you need to specify the visual name for the project and the directory and file to store it.

After you click the *Create* button the [Project View](#) window is opened.