

Captions on Annotations

It is possible to show a value of a qualifier of an annotation instead of the annotation type name in the *Sequence Zoom View*. To enable this option for an annotation type check the *Show value of qualifier* check box and input the values of the required qualifiers in the text field nearby this check box. See the image below.

The screenshot displays a sequence viewer for NC_00136. The main window shows a sequence of amino acids with several annotations. One annotation, a *misc_feature* with a yellow background, is highlighted and contains the text "Show value of qualifier". A red circle highlights this text. Below the sequence is a table with two columns: "Name" and "Value".

Name	Value
▶ CDS	join(2970..3413,3412..3873)
▶ CDS	3875..4999
▶ CDS	5048..5203
▲ misc_feature (0, 2)	
▲ misc_feature	2..590
note	Show value of qualifier
▲ misc_feature	5245..5833
note	Show value of qualifier
▲ source (0, 1)	
▶ source	1..5833

On the right side, the "Annotations Highlighting" panel is visible. It has a section "Select an annotation type:" with a table:

Annotation	Color
CDS	Green
misc_feature	Yellow
source	Grey

Below this is a section "Configure the annotation type:" with three checkboxes:

- Show annotations of this type
- Show on translation
- Show value of qualifier:

The "Show value of qualifier:" checkbox is highlighted with a red circle. Below it is a text input field containing the word "note", also highlighted with a red circle.

If you input several qualifiers names (separated by comma), then the first found qualifier is taken into account and shown on the annotation.